

Product Datasheet

Canine IFNa ELISA Kit (Ready to Use) (orb548706)

Description	Canine Interferon Alpha (IFNa) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for co...
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q75UL5
Sensitivity	5.3pg/mL
Entrez	474181
Expiration Date	Please enquire.