

Product Datasheet

Canine IL2 ELISA Kit (Ready to Use) (orb548709)

Description	Canine Interleukin 2 (IL2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q29416
Sensitivity	6.2pg/mL
Entrez	403989
Expiration Date	Please enquire.