

Product Datasheet

Canine TRAIL ELISA Kit (Ready to Use) (orb548714)

Description	Canine Tumor Necrosis Factor Related Apoptosis Inducing Ligand (TRAIL) Sandwich ELISA KitReady to Us...
Reactivity	Canine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	B4XH66
Sensitivity	0.065ng/mL
Entrez	100174780
Expiration Date	Please enquire.