

Product Datasheet

Gallus IFNa ELISA Kit (Ready to Use) (orb548717)

Description	Chicken Interferon Alpha (IFNa) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for c...
Reactivity	Gallus
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P42165
Sensitivity	28pg/mL
Entrez	100857744
Expiration Date	Please enquire.