

## Product Datasheet

# Gallus IFN beta ELISA Kit (Ready to Use) (orb548718)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Chicken Interferon Beta (IFN beta)<br>Sandwich ELISA Kit<br>Ready to Use: Pre-diluted detection reagents fo...   |
| <b>Reactivity</b>        | Gallus   |
| <b>Range</b>             | 15.6-1000pg/mL   |
| <b>Form/Appearance</b>   | Ready to Use   |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt  |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.  |
| <b>Uniprot ID</b>        | <a href="#">Q90873</a>   |
| <b>Sensitivity</b>       | 5.7pg/mL   |
| <b>Entrez</b>            | <a href="#">554219</a>   |
| <b>Expiration Date</b>   | Please enquire.  |