

Product Datasheet

Other 5-HIAA ELISA Kit (Ready to Use) (orb548725)

Description	General 5-Hydroxyindoleacetic Acid (5-HIAA) Sandwich ELISA Kit Ready to Use: Pre-diluted detection re...
Reactivity	All
Range	6.25-400ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	2.93ng/mL
Expiration Date	Please enquire.