

Product Datasheet

Guinea pig GT ELISA Kit (Ready to Use) (orb548738)

Description	Guinea pig Gastrin (GT) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P06885
Sensitivity	5.9pg/mL
Entrez	100714563
Expiration Date	Please enquire.