

Product Datasheet

Guinea pig INS ELISA Kit (Ready to Use) (orb548740)

Description	Guinea pig Insulin (INS) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Guinea pig
Range	140-9000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P01329
Sensitivity	64.3pg/mL
Entrez	100379579
Expiration Date	Please enquire.