

Product Datasheet

Guinea pig IL2 ELISA Kit (Ready to Use) (orb548745)

Description	Guinea pig Interleukin 2 (IL2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	6.9pg/mL
Expiration Date	Please enquire.