

Product Datasheet

Human FIL1d ELISA Kit (Ready to Use) (orb548825)

Description	Human Interleukin 1 Delta (FIL1d) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9UBH0
Sensitivity	5.7pg/mL
Entrez	26525
Expiration Date	Please enquire.