

Product Datasheet

Human CCL22 ELISA Kit (Ready to Use) (orb548844)

Description	Human Macrophage Derived Chemokine (CCL22) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	000626
Sensitivity	31pg/mL
Entrez	6367
Expiration Date	Please enquire.