

## Product Datasheet

### Human MICB ELISA Kit (orb548847)

<b>Description</b>	Human MHC Class I Polypeptide Related Sequence B (MICB) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Sensitivity</b>	29pg/mL
<b>Expiration Date</b>	Please enquire.