

Product Datasheet

Human Ntn1 ELISA Kit (Ready to Use) (orb548855)

Description	Human Netrin 1 (Ntn1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	095631
Sensitivity	11.6pg/mL
Entrez	9423
Expiration Date	Please enquire.