

Product Datasheet

Human TNFSF13 ELISA Kit (Ready to Use) (orb548894)

Description	Human Tumor Necrosis Factor Ligand Superfamily, Member 13 (TNFSF13) Sandwich ELISA KitReady to Use: ...
Reactivity	Human
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	075888
Sensitivity	16.64pg/mL
Entrez	8741
Expiration Date	Please enquire.