

Product Datasheet

Human VEGFD ELISA Kit (Ready to Use) (orb548906)

Description	Human Vascular Endothelial Growth Factor D (VEGFD) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O43915
Sensitivity	6.1pg/mL
Entrez	2277
Expiration Date	Please enquire.