

Product Datasheet

Monkey IL6 ELISA Kit (Ready to Use) (orb548931)

Description	Monkey Interleukin 6 (IL6) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Monkey
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P79341
Sensitivity	6.4pg/mL
Entrez	102138971
Expiration Date	Please enquire.