

Product Datasheet

Monkey MCP1 ELISA Kit (Ready to Use) (orb548933)

Description	Monkey Monocyte Chemotactic Protein 1 (MCP1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection r...
Reactivity	Monkey
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P61274
Sensitivity	0.054ng/mL
Entrez	102135739
Expiration Date	Please enquire.