

Product Datasheet

Mouse ADIPOR1 ELISA Kit (Ready to Use) (orb548946)

| | |
|--------------------------|--|
| Description | Mouse Adiponectin Receptor 1 (ADIPOR1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent... |
| Reactivity | Mouse |
| Range | 3.12-200ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | Q91VH1 |
| Sensitivity | 1.17ng/mL |
| Entrez | 72674 |
| Expiration Date | Please enquire. |