

Product Datasheet

Mouse ITaC ELISA Kit (Ready to Use) (orb548997)

Description	Mouse Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) Sandwich ELISA Kit Ready to Use: Pre-d...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9JHH5
Sensitivity	12.4pg/mL
Entrez	56066
Expiration Date	Please enquire.