

Product Datasheet

Mouse TGFbR1 ELISA Kit (Ready to Use) (orb549047)

Description	Mouse Transforming Growth Factor Beta Receptor I (TGFbR1) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q64729
Sensitivity	0.069ng/mL
Entrez	21812
Expiration Date	Please enquire.