

Product Datasheet

Mouse TGF beta Receptor 2 ELISA Kit (Ready to Use) (orb549048)

Description	Mouse Transforming Growth Factor Beta Receptor II (TGF beta Receptor 2) Sandwich ELISA KitReady to U...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q62312
Sensitivity	0.105ng/mL
Entrez	21813
Expiration Date	Please enquire.