

Product Datasheet

Mouse TNFSF13 ELISA Kit (Ready to Use) (orb549054)

Description	Mouse Tumor Necrosis Factor Ligand Superfamily, Member 13 (TNFSF13) Sandwich ELISA KitReady to Use: ...
Reactivity	Mouse
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9D777
Sensitivity	22.7pg/mL
Entrez	69583
Expiration Date	Please enquire.