

Product Datasheet

Porcine IGF1 ELISA Kit (Ready to Use) (orb549066)

Description	Porcine Insulin Like Growth Factor 1 (IGF1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection re...
Reactivity	Porcine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P16545
Sensitivity	0.146ng/mL
Entrez	397491
Expiration Date	Please enquire.