

Product Datasheet

Rabbit IGF1 ELISA Kit (Ready to Use) (orb549074)

Description	Rabbit Insulin Like Growth Factor 1 (IGF1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection rea...
Reactivity	Rabbit
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q95222
Sensitivity	0.096ng/mL
Entrez	100008668
Expiration Date	Please enquire.