

## Product Datasheet

### Rat HBEGF ELISA Kit (Ready to Use) (orb549097)

<b>Description</b>	Rat Heparin Binding Epidermal Growth Factor Like Growth Factor (HBEGF) Sandwich ELISA KitReady to Us...
<b>Reactivity</b>	Rat
<b>Range</b>	15.6-1000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q06175</a>
<b>Sensitivity</b>	6.2pg/mL
<b>Entrez</b>	<a href="#">25433</a>
<b>Expiration Date</b>	Please enquire.