

Product Datasheet

Rat CCL22 ELISA Kit (Ready to Use) (orb549107)

Description	Rat Macrophage Derived Chemokine (CCL22) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Rat
Range	23.4-1500pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q5I0L5
Sensitivity	8.45pg/mL
Entrez	117551
Expiration Date	Please enquire.