

Product Datasheet

Bovine EG-VEGF ELISA Kit (Ready to Use) (orb549195)

Description	Bovine Endocrine Gland Derived Vascular Endothelial Growth Factor (EG-VEGF) Sandwich ELISA KitReady ...
Reactivity	Bovine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F1MN36
Sensitivity	6.3pg/mL
Entrez	617599
Expiration Date	Please enquire.