


## Product Datasheet

# Bovine Haptoglobin ELISA Kit (Ready to Use) (orb549212)

<b>Description</b>	Bovine Haptoglobin (Haptoglobin) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for ...
<b>Reactivity</b>	Bovine
<b>Range</b>	156-10000ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q2TBU0</a>
<b>Sensitivity</b>	41.5ng/mL
<b>Entrez</b>	<a href="#">280692</a>
<b>Expiration Date</b>	Please enquire.



Standard Curve of Bovine Haptoglobin ELISA