

Product Datasheet

Bovine HSPB1 ELISA Kit (Ready to Use) (orb549213)

Description	Bovine Heat Shock 27kDa Protein 1 (HSPB1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Bovine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q3T149
Sensitivity	11.9pg/mL
Entrez	516099
Expiration Date	Please enquire.