

Product Datasheet

Bovine HSP90b1 ELISA Kit (Ready to Use) (orb549218)

Description	Bovine Heat Shock Protein 90kDa Beta 1 (HSP90b1) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Bovine
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q95M18
Sensitivity	0.31ng/mL
Entrez	282646
Expiration Date	Please enquire.