

Product Datasheet

Bovine IP10 ELISA Kit (Ready to Use) (orb549236)

Description	Bovine Interferon Gamma Induced Protein 10kDa (IP10) Sandwich ELISA KitReady to Use: Pre-diluted det...
Reactivity	Bovine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q2KIQ8
Sensitivity	0.058ng/mL
Entrez	615107
Expiration Date	Please enquire.