

Product Datasheet

Bovine IL1 alpha ELISA Kit (Ready to Use) (orb549237)

Description	Bovine Interleukin 1 Alpha (IL1 alpha) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Bovine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P08831
Sensitivity	5.8pg/mL
Entrez	281250
Expiration Date	Please enquire.