

Product Datasheet

Bovine LH ELISA Kit (Ready to Use) (orb549253)

Description	Bovine Luteinizing Hormone (LH) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Bovine
Range	0.468-30ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.133ng/mL
Expiration Date	Please enquire.