

## Product Datasheet

# Bovine MIP1a ELISA Kit (Ready to Use) (orb549255)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Bovine Macrophage Inflammatory Protein 1 Alpha (MIP1a) Sandwich ELISA KitReady to Use: Pre-diluted d...   |
| <b>Reactivity</b>        | Bovine  |
| <b>Range</b>             | 15.6-1000pg/mL  |
| <b>Form/Appearance</b>   | Ready to Use  |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.   |
| <b>Uniprot ID</b>        | <a href="#">Q8SQA6</a>  |
| <b>Sensitivity</b>       | 7.4pg/mL  |
| <b>Entrez</b>            | <a href="#">282170</a>  |
| <b>Expiration Date</b>   | Please enquire.   |