

Product Datasheet

Bovine MIP1b ELISA Kit (Ready to Use) (orb549256)

Description	Bovine Macrophage Inflammatory Protein 1 Beta (MIP1b) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Bovine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q17QA1
Sensitivity	9.8pg/mL
Entrez	414347
Expiration Date	Please enquire.