

Product Datasheet

Bovine MCP2 ELISA Kit (Ready to Use) (orb549264)

Description	Bovine Monocyte Chemotactic Protein 2 (MCP2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection r...
Reactivity	Bovine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q09141
Sensitivity	6.4pg/mL
Entrez	281044
Expiration Date	Please enquire.