

## Product Datasheet

# Bovine Neuropeptide Y ELISA Kit (Ready to Use) (orb549267)

<b>Description</b>	Bovine Neuropeptide Y (Neuropeptide Y) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
<b>Reactivity</b>	Bovine
<b>Range</b>	31.2-2000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at - 20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>Q6RUW3</b>
<b>Sensitivity</b>	8.52pg/mL
<b>Entrez</b>	<b>504216</b>
<b>Expiration Date</b>	Please enquire.