

Product Datasheet

Canine ANP ELISA Kit (Ready to Use) (orb549334)

Description	Canine Atrial Natriuretic Peptide (ANP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagen...
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P07499
Sensitivity	6.0pg/mL
Entrez	608289
Expiration Date	Please enquire.