

Product Datasheet

Canine Caveolin 1 ELISA Kit (Ready to Use) (orb549340)

Description	Canine Caveolin 1 (Caveolin 1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for co...
Reactivity	Canine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P33724
Sensitivity	0.064ng/mL
Entrez	403980
Expiration Date	Please enquire.