

Product Datasheet

Canine GMCSF ELISA Kit (Ready to Use) (orb549347)

Description	Canine Colony Stimulating Factor 2, Granulocyte Macrophage (GMCSF) Sandwich ELISA KitReady to Use: P...
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P48749
Sensitivity	7.2pg/mL
Entrez	403923
Expiration Date	Please enquire.