

Product Datasheet

Canine ENG ELISA Kit (Ready to Use) (orb549352)

Description	Canine Endoglin (ENG) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Canine
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	J9P7F0
Sensitivity	0.33ng/mL
Entrez	609166
Expiration Date	Please enquire.