

Product Datasheet

Canine IP10 ELISA Kit (Ready to Use) (orb549366)

Description	Canine Interferon Gamma Induced Protein 10kDa (IP10) Sandwich ELISA KitReady to Use: Pre-diluted det...
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q5KSV9
Sensitivity	6.8pg/mL
Entrez	478432
Expiration Date	Please enquire.