

Product Datasheet

Canine IL17 ELISA Kit (Ready to Use) (orb549372)

Description	Canine Interleukin 17 (IL17) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Canine
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	C6L8D7
Sensitivity	31pg/mL
Entrez	481837
Expiration Date	Please enquire.