

Product Datasheet

Canine LEP ELISA Kit (Ready to Use) (orb549378)

Description	Canine Leptin (LEP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient an...
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	002720
Sensitivity	0.127ng/mL
Entrez	403616
Expiration Date	Please enquire.