

Product Datasheet

Canine LIF ELISA Kit (Ready to Use) (orb549379)

Description	Canine Leukemia Inhibitory Factor (LIF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagen...
Reactivity	Canine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	J9NRL6
Sensitivity	14.7pg/mL
Entrez	403449
Expiration Date	Please enquire.