

## Product Datasheet

# Canine MIP1a ELISA Kit (Ready to Use) (orb549381)

<b>Description</b>	Canine Macrophage Inflammatory Protein 1 Alpha (MIP1a) Sandwich ELISA KitReady to Use: Pre-diluted d...
<b>Reactivity</b>	Canine
<b>Range</b>	78.1-5000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q68A92</a>
<b>Sensitivity</b>	28pg/mL
<b>Entrez</b>	<a href="#">448787</a>
<b>Expiration Date</b>	Please enquire.