

Product Datasheet

Canine MCP1 ELISA Kit (Ready to Use) (orb549386)

| | |
|--------------------------|---|
| Description | Canine Monocyte Chemotactic Protein 1 (MCP1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection r... |
| Reactivity | Canine |
| Range | 0.156-10ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | P52203 |
| Sensitivity | 0.058ng/mL |
| Entrez | 403981 |
| Expiration Date | Please enquire. |