

Product Datasheet

Canine NT-ProBNP ELISA Kit (Ready to Use) (orb549389)

| | |
|--------------------------|--|
| Description | Canine N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) Sandwich ELISA Kit Ready to Use: Pre-dilu... |
| Reactivity | Canine |
| Range | 78.1-5000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Sensitivity | 28pg/mL |
| Expiration Date | Please enquire. |