

Product Datasheet

Canine iNOS ELISA Kit (Ready to Use) (orb549390)

Description	Canine Nitric Oxide Synthase 2, Inducible (iNOS) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O62699
Sensitivity	0.126ng/mL
Entrez	403822
Expiration Date	Please enquire.