

Product Datasheet

Canine TSH ELISA Kit (Ready to Use) (orb549401)

Description	Canine Thyroid Stimulating Hormone (TSH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Canine
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	11.25pg/mL
Expiration Date	Please enquire.